

Procedure for Immunofluorescence (IF) Staining of Slides with Enriched Populations of Circulating Tumour Cells (CTCs)

1. MATERIALS, EQUIPMENT AND FORM

| Reagents & Solutions: | Equipment: | Supplies: |
|---|---|------------------------------|
| Skim milk powder | Humidified chamber | Coverslips (22x40 and 22x50) |
| Tween-20 | Micropipettors | 10, 100 and 1000 µl tips |
| 1x PBS | Coplin jars | Laboratory wipes |
| Vectashield Mounting Medium with DAPI | Magnetic stirrer | |
| CD45 Primary antibody: AbCam Rabbit Monoclonal [EP322Y] to CD45 (Cat. # ab40763) | Stir bar | |
| Texas Red Secondary antibody: Sheep Polyclonal Secondary Antibody to Rabbit IgG – H&L [Texas Red] (Cat. # ab6793) | Aqueous waste container | |
| Pan-CK antibody: WaveSense Pan Cytokeratin – FITC (4,5,6,8,11,13,18) Cat. # R2105-1 | Fluorescence microscope equipped with recommended filters | |

2. REAGENTS

1x PBS

Add 100mL of 10x PBS to 900mL of distilled water.

Wash buffer

Add 1mL Tween-20 to 999mL 1x PBS, mix well.

Prepare fresh (can be stored for 2 weeks at room temperature).

5% Skim milk

Add 5g skim milk powder to 100mL wash buffer.

Place on a stirrer until the skim milk dissolves.

Always make fresh.

3. PROCEDURES

SAMPLES

- Samples consist of patient blood samples, or donor blood/buffer samples spiked with cancer cells from culture, and applied to slides using the Wavesense immunomagnetic enrichment assay.
- Following fixation, slides should be kept in a desiccated slide box at -20°C until ready for the IF assay.



QUALITY CONTROL

- Control slides must be run concurrently with patient slides to monitor assay performance and to assess the accuracy of signal enumeration.
- Cell line or normal blood spreads should be used as positive and negative controls for the IF testing.
- Controls should be run on each day of IF testing.

PREPARATIONS REQUIRED BEFORE STARTING PROCEDURE

- All reagents and stock solutions should be prepared prior to the start of the procedure
- Label the slides correctly: patient ID, date, lab identification, etc., along with name of person doing the procedure
- Prepare fresh solutions prior to each procedure
- To protect the light sensitive fluorophores, work in reduced light conditions as much as possible

PROTOCOL

Treating slide with blocking agent

- Remove the slides from -20°C freezer and bring them to room temperature.
- Rinse the slide with 1x PBS briefly in a coplin jar.
- Place the slide in a coplin jar containing 5% skim milk and incubate at room temperature for 1 h.
- Wash the slide 5 times with RT 1x PBS/0.1% Tween-20.
- Hold the slide vertically on a kimwipe in order to drain entire wash buffer, do not allow slide to dry.

Primary CD45 antibody assay

- Make a 1:200 dilution of Rabbit and Human CD45 primary antibody in Dako antibody diluent.
- Apply 100µl of primary antibody on to the slide.
- Cover the slide with a 22x50mm coverslip and place in a humidified chamber at room temperature for 45 min.
- Gently remove the coverslip and wash the slide 5 times with RT 1x PBS/0.1% Tween-20 in a coplin jar replacing fresh wash buffer after every wash.
- Hold the slide vertically on a kimwipe in order to drain the entire wash buffer, do not allow slide to dry.

Secondary CD45 antibody assay

- Make a working solution of 2.5 µg/mL Sheep anti-rabbit polyclonal antibody (Texas Red) in Dako antibody diluent.
- Apply 100ul of secondary antibody on to the slide.
- Cover the slide with a 22x50mm coverslip and place in a humidified chamber at room temperature for 15 min in dark.
- Gently remove the coverslip and wash the slide 5 times with RT 1x PBS/0.1% Tween-20 in a coplin jar replacing fresh wash buffer after every wash.
- Hold the slide vertically on a kimwipe in order to drain the entire wash buffer, do not allow slide to dry.

PanCK antibody assay

- Apply 100µl Mouse anti-human PanCK FITC conjugate to the slide.
- Cover the slide with a 22x50mm coverslip and place in a humidified chamber at room temperature for 15 min in dark.
- Gently remove the coverslip and wash the slide 5 times with RT 1x PBS/0.1% Tween-20 in a coplin jar replacing fresh wash buffer after every wash.
- Hold the slide vertically on a kimwipe in order to drain the entire wash buffer, do not allow slide to dry.

Counterstain assay

- Apply 20µl of Vectashield Mounting Medium with DAPI to the target area and apply coverslip (22x50)
- Keep the slides in a desiccated slide box at -20°C for at least 30 minutes.
- View the slides under appropriate filters under a fluorescence microscope.

Table 2. Summary Protocol

| | Action | Temperature | Incubation time |
|----|--|--------------------|------------------------|
| 1 | Remove the slide from -20 C freezer and bring them to room temperature | RT | |
| 2 | Rinse the slide with 1x PBS briefly in a coplin jar | RT | |
| 3 | Incubate the slide in a coplin jar containing 5% skim milk | RT | 60 minutes |
| 4 | Wash the slide 5 times with 1x PBS/0.1% Tween-20 | RT | 2 minutes each |
| 5 | Apply 100 µl of primary antibody on to the slide | RT | |
| 6 | Cover the slide with a 22x50mm coverslip | RT | |
| 7 | Incubate slide in a humidified chamber | RT | 45 minutes |
| 8 | Wash the slide 5 times with 1x PBS/0.1% Tween-20 | RT | 2 minutes each |
| 9 | Apply 100 ul of secondary antibody on to the slide | RT | 60 minutes |
| 10 | Cover the slide with a 22x50mm coverslip | RT | |
| 11 | Incubate slide in a humidified chamber in the dark | RT | 15 minutes |
| 12 | Wash the slide 5 times with 1x PBS/0.1% Tween-20 | RT | 2 minutes each |
| 13 | Apply 100µl Mouse anti-human PanCK FITC conjugate to the slide | RT | |
| 14 | Cover the slide with a 22x50mm coverslip | RT | |
| 15 | Incubate slide in a humidified chamber in the dark | RT | 15 minutes |
| 16 | Wash the slide 5 times with 1x PBS/0.1% Tween-20 | RT | 2 minutes each |
| 17 | Apply 20 µl Vectashield Mounting Medium with DAPI to the target area and apply coverslip (22x50mm) | RT | |
| 18 | Keep the slides in a desiccated slide box at -20°C Freezer | RT | 30 minutes |
| 19 | View the slides using appropriate filters under a fluorescence microscope | | |

References

- Refer to the 3rd edition of the Public Health Agency of Canada's Laboratory Biosafety Guidelines,



- Refer to the Queen's University's Biohazards Safety Manual published by the Department of Environmental Health and Safety, when handling biohazardous materials such as blood and tumor samples