

Directions for AQIX® Cryopreservation of HS-cells [SOP/HS/CRYO]

- Step 1.** Harvest HS-cells at 10^7 to 10^6 cells/mL.
- Step 2.** Prepare pellets of HS-cells at 10^7 to 10^6 cells/mL by gentle centrifugation in suitable vials.
- Step 3.** Place HS-cell pellet vials, unopened containers of sterile AQIX® RS-I/1G solution and cryopreservative vials (e.g., NUNC™ Cryoflex) into a laminar flow cabinet.
- Step 4.** *Gently* add the cryopreservative fluid, AQIX® RS-I/1G to the HS-cell pellets to achieve 10^7 to 10^6 cells/mL suspension volumes of the cryopreserve fluid to a maximum of 1-2 mL volumes.
- Step 5.** Allow HS-cell suspensions to equilibrate at ambient temperatures for 15 – 30 minutes within the laminar flow cabinet.
- Step 6.** Place 1-2 mL volumes of the HS-cell suspension media into suitable cryopreservation vials (e.g., NUNC™ Cryoflex) and seal tightly to ensure that no leakage of Nitrogen gas into the vials occurs during cooling and storage.
- Step 7.** Remove sealed, cryopreservation vials from the laminar flow cabinet and transfer to a suitable cryopreservation system (e.g., Nalgene “Mr. Frosty”; Cat.No. 5100-0001).
- Step 8.** Cool HS-cell suspension vials at $1^\circ\text{C}/\text{min}$ to -40°C and then transfer into a suitable storage unit to achieve the final, desired cryopreservation temperature (e.g., -60°C to -150°C).
- Step 9.** After 48 hours of cryopreservation, remove a sample vial to test for cell viability. Proceed as directed in Step 10.
- Step 10.** To reconstitute cryopreserved HS-cells, allow the frozen vials of suspension media to defrost at 37°C in a tissue culture cabinet or (less preferred) within a conventional water bath.
- Step 11.** *Gently* pipette off HS-cell suspension and transfer 1.0 mL of the suspension into either (a) 2.0 mL of (1x) AQIX® RS-I solution or (b) 2.0 mL of a preferred TC-media. Repeat this procedure a further two times to dilute out the AQIX® RS-I/1G cryopreservation fluid.
- Step 12.** Test HS-cell viability using conventional procedures.