

[SOP/H-UC/ST/11-08]

## ***Human Umbilical Cord Collection Protocols:***

### **Method:**

#### **Step 1. – Umbilical Cord Procurement Site Preparation**

Prepare a clean, sterile surface area in the ante-room to the operating theatre.

Remove the bottle of AQIX<sup>®</sup> UC fluid and swab the outside of the bottle with a suitable sterilising agent. Assemble the bottle of AQIX<sup>®</sup> UC fluid along with surgical gloves, face mask, scissors and forceps on the prepared surface.

#### **Step 2. – Collecting the Umbilical Cord**

##### **Sample Retrieval Procedures:**

- a) Put on surgical gloves and face mask and retrieve the placenta plus umbilical cord biopsy from operating theatre and transfer to the prepared area selected.
- b) Using the sharp scissors, carefully cut the cord about 5-10cm from the placenta.
- c) Holding the umbilical cord (UC), cut through the remaining placenta leaving a 5 - 10cm circular rim of the outer, epithelial (amnion) layer containing multipotent cell types (Miki et al, 2005).
- d) Using forceps, mechanically peel off the chorion from the amnion layer to form the UC-amnion sample.
- e) Thoroughly rinse the 'UC- amnion' biopsy several times using AQIX<sup>®</sup> UC-fluid and finally store/transport the UC-amnion HSC sample in a 250 mL bottle of AQIX<sup>®</sup> RS-I **or** its antibiotic (RS-I/Ab) **or** antifungal (RS-I/Af) **or** combination (RS-I/Ab+Af) formulations.
- f) Screw lid firmly onto the UC-sample bottle and ensure that there are no leakages.
- g) Place AQIX<sup>®</sup> UC-sample bottle into Transport Box for dispatch to a laboratory for further processing.

#### **Step 3 – Storage & Transportation Procedures:**

##### **A) Hypothermic Conditions (recommended)**

- l) Place a **frozen** cool pack on the bottom of the polystyrene transportation box and the other '**cold**' but unfrozen cool pack on top of the solid, **frozen** cool pack to act as a 'cushion' for the AQIX<sup>®</sup> UC-sample bottle.

- m) Place the AQIX<sup>®</sup> UC-sample bottle on top of the *unfrozen* cool pack and fill the surrounding spaces with sterile packaging to stop the sample bottle moving about and to retain a low temperature during transport.
- n) Place the lid on the transportation container and tape closed the cardboard sleeve to ensure contents are safe.
- o) Dispatch transportation to recipient laboratory for isolation of placenta (amnion) and chord HSCs.

**A) Ambient Conditions**

- l) Place the AQIX<sup>®</sup> UC-sample bottle into a polystyrene transportation box and fill the surrounding spaces with sterile packaging to stop the sample bottle moving about and to retain a low temperature during transport.
- n) Place the lid on the transportation container and tape closed the cardboard sleeve to ensure contents are safe.
- o) Dispatch transportation to recipient laboratory for isolation of placenta (amnion) and chord HSCs.

**Step 4 – Isolation of Umbilical Cord (UC) and Amnion Epithelium (AE) Derived Cells:**

- p) Remove the sample from transportation container and;
  1. Either, Prepare UC-tissue for immediate extraction and culture of stem cells as described by Miki et al., 2005 [Stem Cells.2005;23(10):1549-59]
  - Or,
  - 2 Subdivide UC-tissue for short-term preservation (1 - 5 days) in 'fresh' AQIX RS-I @ 3-8 °C
  - Or,
  - 3 Freeze gradually over 16-24 hours @ **-60** °C for Long-term storage using a suitable cryopreservation solution (DMSO; AQIX RS-I + 5% Glycerol)

Author: Dr. D. Rees MRSNZ  
CSO  
Aqix Ltd (UK)