

## Blood Processing

### A. Reagents Needed

DNAzol  
PBS Sterile solution 1X  
Ficoll  
ACK Lysing Buffer  
Trypan Blue  
Chilled Filtered 0.22 $\mu$ m 100% FCS  
Chilled Filtered 0.22 $\mu$ m 80% FCS 20% DMSO

**B. Objective:** to separate and store DNA, Blood Plasma and PBMCs.

### C. Description

#### DNA

1. Invert the blood tube 2-3 times to ensure homogenous mixture.
2. Using a sterile pipet, take out 1ml of whole blood and transfer into 15ml falcon tube.
3. Mix whole blood with 2ml of DNAzol.
4. Aliquot Blood-DNAzol mixture (1.5ml/vial) into two labeled cryovials.
5. Store vials immediately at -80° C.

#### Plasma

1. Centrifuge blood samples at 1500RPM for 10 minutes at Room temperature (RT).
2. Reserve spun blood.
3. Harvest supernatant (Plasma) into separate 15ml falcon tube.
4. Re-spin the plasma at 2500RPM for 5 minutes at RT.
5. Aliquot supernatant into 5 labeled round-bottom cryovials and store in -80° C

#### Isolating PBMC

6. Add 5ml of Ficoll to each 15ml tube (or 10 of Ficoll to the 50ml Falcon tube).
7. Re-suspend the whole blood to two times its original volume with sterile PBS.
8. Homogenize by pipetting up and down.
9. Use a 15ml Falcon tube if the amount of blood is below 4.5ml, if above use 50ml tubes.

10. Tip tube of Ficoll almost horizontally and add diluted blood **EXTREMELY SLOWLY** right above the Ficoll. Do not let the blood and Ficoll mix, Ficoll and blood layers should remain separate.
11. Centrifuge the Ficoll/blood at 2500RPM for 15 minutes at RT. **IMPORTANT: set deceleration to "0" and acceleration to "3"**
12. Harvest the cell ring (PBMC) located above the Ficoll in a circular movement.
13. Transfer cells into a second 15ml or 50ml tube filled with PBS (5ml for 15ml tube/15ml for 50ml tube).
14. Spin cells at 1500RPM for 5 minutes (spin 10min if using 50ml tube) at RT.
15. Discard supernatant and gently tap bottom of tube to re-suspend pellet.
16. Re-suspend cells again with PBS 1X (if pellet looks red, re-suspend cells in 3-5ml of ACK buffer to lyse remaining RBC's, wash with PBS after 5 minutes) Transfer the cell suspension into a 15ml tube, if originally in 50ml tube.
17. Spin cells again at 1500RPM for 5 min at RT.
18. Discard supernatant, re-suspend in 5ml PBS.

### Cell Counting

19. Pipet cell suspension up and down to homogenize. Take a drop from the center.
20. Pipette 10 $\mu$ l of Trypan Blue into a test tube, add 10 $\mu$ l of cell suspension and mix thoroughly. Add 10 $\mu$ l of cell/dye mixture to hemocytometer.
21. Count 2 quadrants of cells, take the average, and apply cell concentration formula: #of cells x 2 x total volume of cell suspension x 10<sup>4</sup> (e.g. 40 cells counted in 5ml of PBS  $\rightarrow$  40 x 2 x 5 x 10<sup>4</sup> = 4 x 10<sup>6</sup> cells total). Count at least 50 cells, if less than 50 cells in 2 quadrants, count all 4 quadrants and average.

### Cryopreservation

22. Determine number of round-bottom cryovials to make as follows:
  - Less than 3 x 10<sup>6</sup> cells: 2 cryovials.
  - 3-4 x 10<sup>6</sup> cells: 3 cryovials.
  - 4-5 x 10<sup>6</sup> cells: 4 cryovials.
  - Continue up to 5 cryovials max post-transplant; 15 max pre-transplant (amount may vary depending on study/organ type).
23. FREEZE ONLY 10 or less aliquots at a time.
24. Always keep 0.22 $\mu$ m filtered FCS and FCS/DMSO in ice.
25. After counting cells, centrifuge cell suspension at 1500RPM for 5 minutes at RT.
26. Label appropriate number of cryovials and chill on ice.
27. Discard supernatant and re-suspend cells in 0.22 $\mu$ m filtered FCS (0.5ml/vial to be made).

28. Distribute 0.5ml of cell re-suspended in ice-cold FCS into each chilled cryovial.
29. Add 0.5ml ice cold 0.22 $\mu$ m filtered 80% FCS/ 20% DMSO to each cryovial of cells to make a 1:1 solution of FCS/FCS-DMSO
30. Recap and turn vial containing cells and preservative upside down several times to homogenize and quickly store in -80° C.
31. Long-term storage of frozen cells in Li N<sub>2</sub>.