

Date: September 6, 2022

Department of Translational Molecular Pathology (TMP)

Immune Profiling Laboratory

Cancer Immune Monitoring and Analysis Center


The University of Texas MD Anderson Cancer Center

TMP-Immunoprofiling Laboratory Director:

Cara Haymaker, Assistant Professor, TMP

Immunohistochemistry and Digital Pathology Lab Director:

Luisa Maren Solis Soto, TMP



SOP for Pathological Assessment and Region of Interest Selection using the nanoString GeoMx Digital Spatial Profiler

1. Technical platform

The GeoMx® Digital Spatial Profiling (DSP) technology is a novel platform developed by NanoString. This product relies upon antibodies or nucleic acid probes coupled to photocleavable oligonucleotide tags. After the incubation of antibodies to slide-mounted tissue sections, cell types of interest with fluorescent morphology markers to create fluorescent images for selecting region of interest (ROIs). The oligonucleotide tags are released from discrete regions of the tissue via UV exposure. Released tags are quantitated in an nCounter® and counts are mapped back to tissue location, yielding a spatially resolved digital profile of analyte abundance(1).

This SOP describes the **workflow and process for experimental design, pathological assessment and region of interest selection** using the GeoMx DSP(2) at the Translational Molecular Pathology Immunoprofiling Laboratory (TMP-IL).

2. Experimental design

Profiling tissue using the GeoMx Digital Spatial profiler requires several considerations to design and conduct DSP assay in order to reach scientific aims of a specific project. To tailor selection of visualization biomarkers, and selection of DSP panels and modules, the following information should be obtained during scientific discussion with collaborators. (See also **Annex 1**. TMP Experimental design form)

a. Sample information:

- i. Number of samples: Number of FFPE blocks
- ii. Number of slides: Specify if multiple samples can be placed in one slide or samples can be placed in TMA.
- iii. Tumor type: Histopathology diagnosis of samples.
- iv. Type of samples: Surgically resected, core biopsies, TMA, decalcified samples or others.

b. Aims:

- i. **Principal aim.**
- ii. **Specific Scientific aims.**

c. GeoMx DSP biomarker/panel and ROI design:

- i. **Morphology biomarkers:** Determine if the project will use commercially available GeoMx DSP morphology kits or if the laboratory will need to optimize a custom morphology biomarker

panel. A validation process will be performed for custom panels. The following morphology biomarkers are available and optimized at TMP-IL:

| Markers | Clone | Vendor | Catalog # | Labeling Fluorophore |
|---------------|-------------|--------------------|-----------------|----------------------|
| PanCK | AE1+AE3 | Nanostring/Novus | NBP2-33200 | AF532 |
| CD45 | 2B11+PD7/26 | Nanostring/CST | #13917 | AF594 |
| S100B | S100B/1706R | Nanostring / Novus | NBP2-54426 | AF532 |
| B3 Tubulin | EP1569Y | abcam | ab190575 | AF647 |
| CD20 | IGEL/773 | NovusBio | NBP2-47840 | AF594 |
| CD3e | UMAB54 | Origene | AC211030 | AF647 |
| CD68 | KP-1 | Santa Cruz | SC-20060 AF647 | AF647 |
| CD68 | KP-1 | Santa Cruz | sc-20060 AF594 | AF594 |
| Claudin-4 | 3E2C1 | LSBio | LS-C354893 | AF594 |
| NFH | EPR20020 | abcam | ab223139 | AF594 |
| PanCK | AE1/AE3 | NovusBio | NBP2-33200AF647 | AF647 |
| Synaptophysin | YE269 | abcam | ab 196166 | AF647 |

- ii. **DSP profiling panel:** Protein or RNA.
- iii. **DSP core panel and modules:** The following core panels and modules are currently commercially available for human samples (See also **Annex 2**).

Protein:

- **Immune Cell Profiling Core (Human protein core):** Includes key immuno-oncology targets and markers of immune cell types, including T cells, B cells, macrophages, NK cells, epithelia, and stroma. It also contains the controls needed to run any GeoMx DSP experiment.

Biomarkers: Beta-2-microglobulin; CD3; CD56; CTLA4; GZMB; PD-1; CD11c; CD4; CD68; Pan-cytokeratin; HLA-DR; PD-L1; CD20; CD45; CD8; Fibronectin; Ki-67; SMA; Ms IgG2a; Ms IgG1; Rb IgG; S6; Histone H3; GAPDH.

- **IO Drug Target Module:** Includes drug targets in development within the immuno-oncology space, including checkpoint molecules and metabolic mediators of immune function.

Biomarkers: 4-1BB; LAG3; ARG1; OX40L; B7-H3; STING; GITR; TIM-3; IDO1; VISTA

- **Immune Activation Status Module:** Includes additional checkpoint molecules that modulate T cell activation.

Biomarkers: CD127; CD25; CD27; CD40; CD44; CD80; ICOS; PD-L2

- **Immune Cell Typing Module:** Includes an expanded set of cell type markers to more deeply profile immune cell types covered in the Immune Cell Profiling Core and measure additional immune cell types, including T cell subsets.

Biomarkers: CD14; CD163; CD34; CD45RO; CD66b; FAPalpha; FOXP3

- **Pan-Tumor Module:** Includes an expanded set of markers for detecting specific tumor types, including ER+/HER2+ breast tumors, hematopoietic malignancies, and melanoma.

Biomarkers: BCL-2; EpCAM; ER alpha; HER2/ERBB2; MART1; NY-ESO-1; PR; PTEN; S100B

• *Cell Death Module:* Includes protein mediators of immunogenic and programmed cell death.
Biomarkers: BAD; BCL6; BCLXL; BIM; CD95/Fas; GZMA; NF1; p53; PARP; Cleaved Caspase 9

• *PI3K/AKT Signaling Module:* Includes key proteins involved in PI3K-AKT signal transduction, and phosphorylated protein products that measure pathway activation.

Biomarkers: Pan-AKT; MET; Phospho-AKT1(S473); Phospho-GSK3B (S9); Phospho-Tuberin (T1462); Phospho-GSK3A (S21)/Phospho-GSK3B (S9); INPP4B; PLCG1; Phospho-PRAS40 (T246)

• *MAPK Signaling Module:* Includes key proteins involved in MAPK signal transduction, and phosphorylated protein products that measure pathway activation

Biomarkers: EGFR; pan-RAS; BRAF; Phospho-c-RAF (S338); Phospho-JNK (T183/Y185); Phospho-MEK1 (S217/S221); Phospho-p38 MAPK (T180/Y182); Phospho-p44/42 MAPK ERK1/2 (T202/Y204); p44/42 MAPK E RK1/2; Phospho-p90 RSK (T359/S363)

RNA:

- GeoMx Immune pathways panel: 84-plex RNA panel + customs (96-plex)/nCounter
- Cancer Transcriptomic Atlas (CTA) / Next Generation Sequencing: Covers +1,800 genes
- Whole Transcriptomic Atlas (WTA)/ Next Generation Sequencing: Covers +18,000 genes

d. Region of Interest (ROI) strategy

- i. Number of ROI per sample.
 - ii. Shape of ROI: rectangles, circles and/or polygons up to 660x785 μm .
 - iii. Histological or phenotypic strategy for ROI placement.
 - iv. Segmentation of ROI: The strategy should state which morphology biomarkers will be used to segment the area.
- e. **Estimated time to complete the project.** For Protein and 84plex RNA panels the maximum capacity is 8 slides per week at 12 AOIs per sample (Total 96 AOIs per week)

3. Selection of cases and tissue qualification for DSP analysis

After a cohort of cases has been selected for DSP analysis, a pathological assessment of samples available is performed using freshly cut sections stained with H&E (preferable) or archived H&E slides, the latter is less preferable due to the risk of evaluation of a section that is not representative of tissue left in FFPE block. (See also Annex 3. check list for selection of cases and tissue qualification for DSP analysis)

- a. **Pathology assessment of samples:** The laboratory requires that the H&E slides of potential tissue samples for DSP assay are evaluated by a pathologist prior to sending material for DSP assay. The H&E slides are digitally scanned using Aperio AT2 (Leica Biosystems). An excel file from collaborators containing relevant clinicopathological information of the samples should also be provided.

The following information is usually provided by collaborators but may vary among projects.

- MRN
- Tissue ID
- Time Point. If the patient has two TPs (TP0, TP1), check the experimental design for specific

instructions on sample selection.

- Diagnosis
- Organ site.
- Type of biopsy (surgical specimen, core biopsy, endoscopic biopsy, TMA) (A map of the TMA must be attached)
- Localization (primary organ vs metastasis in lymph nodes or other organs).
- Other

b. **Quality control:** A pathologist will perform a histological quality control of samples received using H&E slides (scanned and/or physical slide). The following information is usually evaluated for surgically resected tumor samples and can be tailored for specific projects. The following information is evaluated in surgical resected tumor samples

- Presence of tumor cells (Yes/No): > than 100 viable tumor cells.
- Presence of normal adjacent tissue.
- Presence of Tumor associated Tertiary lymphoid aggregates (TLS).
- Presence of Invasive margin (IM).
- Percentage of tumor infiltrating lymphocytes (TILs)*.
- Comments: Annotations of potential presence of artifacts that can interfere with the staining (background) including blood cells, fibrin, bone or necrosis.

*Note: TIL's assessment is performed in each sample to assess % of immune cells in tumor stroma using the guidelines of the International TILs Working Group 2014 (4).

In samples with heterogeneous tumor immune infiltration, we use a grid of 1500 μm^2 to evaluate TILs% in each square to obtain an average of TILs %; this evaluation can be further used to select ROI for DSP analysis.

c. **Tissue qualification for DSP assay:** After histological evaluation a pathologist will perform the following:

- i. Select the best representative sample to run DSP assay (Considering experimental design strategies).
- ii. Annotate potential ROI areas (squares, circles or polygons of up to $660 \times 785 \mu\text{m}^2$) according to the experimental design strategy.

d. **Requisition of material to perform DSP assay:** After a pathologist have provided a list with the selected samples, the laboratory will request freshly-cut unstained sections processed according to nanoString GeoMx DSP protocol(5) to be sent to TMP-IL along with excel file containing information of slides.

e. **Storage of slides:** Once the laboratory receives unstained sections and a file with sample identification.

- i. Sections will be stored at 4°C.
- ii. Slides must be labeled with sample identification that match with information of slides submitted in the excel file.
- iii. Slides labels should not contain protected health information (PHI).
- iv. Date of sectioning should be provided in the slide label, slide container or in an insert inside the container.

- f. **Quality assurance of unstained sections:** Unstained sections are reviewed by a pathologist to assure that:
- i. The sections are adequate with no excessive folding or detachment
 - ii. Tissue fit the DSP scanning area
 - iii. No tissue overlap the gasket area or the tip calibration area.
 - iv. If multiple unstained slides are available to select for this assay, a pathologist will choose the unstained section with best tissue integrity.
 - v. Unstained sections should be similar in shape to H&E previously evaluated. If significant differences are found, H&E of a serial sections should be requested.

4. Selection of ROI in the scanned images using DSP:

Once the slides have being processed and the slides have been scanned in the DSP device, a pathologist will access the GeoMx DSP software through Google Chrome based app, via MD Anderson network, and will login using the username and password previously set up by the laboratory. Using the GeoMx DSP software tools, a pathologist will perform the following tasks (See also **Annex 4** Selection of Regions of interest and segmentation using GeoMx DSP software):

- i. Visualize the slide image which will allow to perform a verification of the immunofluorescence staining quality of morphology biomarkers, and
- ii. Select ROI and perform segmentation according to experimental design strategy.
- iii. Approve ROI and segment and send for collection.

5. References

1. Van TM, Blank CU. A user's perspective on GeoMx™ digital spatial profiling. *Immuno-Oncology Technology*. 2019;1:11-8.
2. Hernandez S, Lazcano R, Serrano A, Powell S, Kostousov L, Mehta J, et al. Challenges and Opportunities for Immunoprofiling Using a Spatial High-Plex Technology: The NanoString GeoMx(®) Digital Spatial Profiler. *Front Oncol*. 2022;12:890410.
3. NanoString I. GeoMx® Protein Assays for Immuno-Oncology 2021 [cited 2022 09/05/2022]. Available from: chrome-extension://efaidnbmnribpcajpcglclefindmkaj/https://nanostring.com/wp-content/uploads/PB_MK3350_DSP_IO_Protein_PB_R28.pdf.
4. Salgado R, Denkert C, Demaria S, Sirtaine N, Klauschen F, Pruneri G, et al. The evaluation of tumor-infiltrating lymphocytes (TILs) in breast cancer: recommendations by an International TILs Working Group 2014. *Ann Oncol*. 2015;26(2):259-71.
5. NanoString I. GeoMx® DSP Automated Slide Preparation User Manual 2022 [updated April 2022; cited 2022 06/27/2022]. Available from: <https://university.nanostring.com/geomx-dsp-automated-slide-preparation-user-manual/1209595>.

Immunohistochemistry and Digital pathology Laboratory

GeoMx DSP scientist leads:

Sharia Hernandez, MD, Research scientist.

Alejandra Serrano, Research scientist.

Director: Luisa Maren Solis Soto, MD Assistant Professor MD Anderson Cancer Center

Signature (Luisa M. Solis Soto)

Date: 09/06/2022

Annex 1. GeoMx DSP experimental design form.

**GEOMX DSP EXPERIMENTAL DESIGN FORM
DSP project # (n= number of slides)**

- 1) **Name of the project:**
 - 2) **Project PI:**
 - 3) **Cohort:**
 - 4) **Sample information:**
 - i) Number of samples: *Number of FFPE blocks*
 - ii) Number of slides: *Specify if multiple samples can be placed in one slide or samples can be placed in TMA.*
 - iii) Tumor type: *Histopathology diagnosis of samples.*
 - iv) Type of samples: *Surgically resected, core biopsies, TMA, decalcified samples or others.*
 - b) **Aims**
 - i) **Principal aim:**
 - ii) **Specific Scientific aims:**
 - c) **GeoMx DSP biomarker/panel and ROI design:**
 - i) **Morphology biomarkers:** *Determine if the project will use commercially available GeoMx DSP morphology kits or if the laboratory will need to optimize a custom morphology biomarker panel. A validation process will be performed for custom panels.*
 - ii) **DSP profiling panel:** *Protein or RNA.*
 - iii) **DSP core panel and modules:** *The following core panels and modules are currently commercially available: (See list in Annex 2)*
 - d) **Region of Interest (ROI) strategy**
 - i) Number of ROI per sample:
 - ii) Shape of ROI: *rectangles, circles and/or polygons up to 660 x 785 μ m.*
 - iii) Histological or phenotypic strategy for ROI placement:
 - iv) Segmentation of ROI: *The strategy should state which morphology biomarkers will be used to segment the area, order of collection and segment definitions*
 - i. Segment 1:
 - (a) Segment definition:
 - ii. Segment 2:
 - (a) Segment definition
 - iii. Segment 3:
 - (a) Segment definition
 - e) **ROI strategy for tonsil control if included:**
 - i) Number of ROI: 9 ROIs.
 - ii) Shape of ROI: polygons up to 660 x 785 μ m
 - iii) Histological or phenotypic strategy for ROI placement: 3 ROIs will be placed in germinal center; 3 ROI in interfollicular area; 3 ROI in reticulated epithelium area. No segmentation.
 - f) **Plan to select ROIs per week:**
 - i) 1st week: 8 cases
 - ii) 2nd week: 7 cases + control-Tonsil
 - iii) Estimated number of ROIs:
 - iv) Estimated number of AOIs:
- 5) **TMP-IL Personnel**

SOP for Pathological Assessment and Region of Interest Selection using the nanoString GeoMx DSP

- a) GeoMx DSP Platform Operator:
 - b) GeoMx DSP Platform Lab member:
 - c) Sectioning:
 - d) GeoMx DSP Pathologists:
 - e) TMP-IL IHC /DSP Lab director:
- 6) **Estimated time to complete the project.** For Protein and 84plex RNA panels the maximum capacity is 8 slides per week at 12 AOI per sample (Total 96 AOI per week)

-----END-----

Annex 2 GeoMX DSP panels commercially available for human specimens as for 09/06/2022

1. PROTEIN

| Nanostring DSP Protein Immuno-oncology Panels | | | |
|--|--|--------------------------|--|
| Immune Cell Profiling | IO Drug Target Module | Cell Death Module | PI3K/AKT Signaling Module |
| Core | | | |
| PD-1 | 4-1BB | BAD | Pan-AKT |
| CD68 | LAG3 | BCL6 | MET |
| HLA-DR | OX40L | BCLXL | Phospho-AKT1(S473) |
| Ki-67 | Tim-3 | BIM | Phospho-GSK3B (S9) |
| Beta-2-microglobulin | VISTA | CD95/Fas | Phospho-Tuberin (T1462) |
| CD11c | ARG1 | GZMA | Phospho-GSK3A (S21)/Phospho-GSK3B (S9) |
| CD20 | B7-H3 | NF1 | INPP4B |
| CD3 | IDO1 | p53 | PLCG1 |
| CD4 | STING | PARP | Phospho-PRAS40 (T246) |
| CD45 | GITR | Cleaved Caspase 9 | |
| CD56 | | | MAPK Signaling Module |
| | Immune Activation Status Module | | EGFR |
| CD8 | CD127 | | pan-RAS |
| CTLA4 | CD25 | | BRAF |
| GZMB | CD80 | | Phospho-c-RAF (S338) |
| PD-L1 | ICOS | | Phospho-JNK (T183/Y185) |
| PanCk | PD-L2 | | Phospho-MEK1 (S217/S221) |
| SMA | CD40 | | Phospho-p38 MAPK (T180/Y182) |
| Fibronectin | CD44 | | Phospho-p44/42 MAPK ERK1/2 (T202/Y204) |
| Rb IgG | CD27 | | p44/42 MAPK E RK1/2 |
| Ms IgG1 | | | Phospho-p90 RSK (T359/S363) |
| Ms IgG2a | | | |
| | Immune Cell Typing Module | | |
| Histone H3 | CD45RO | | |
| S6 | FOXP3 | | |
| GAPDH | CD34 | | |
| | CD66b | | |
| | FAP-alpha | | |
| | CD14 | | |
| | CD163 | | |

2. RNA:

2.1. GeoMx Immune Pathway Panel Human RNA Core for nCounter (84-plex+custom, 96-plex)

| Probe Name | Gene(s) |
|-------------------|----------------|
| AKT1 | AKT1 |
| ARG1 | ARG1 |
| B2M | B2M |
| BATF3 | BATF3 |
| BCL2 | BCL2 |
| CCL5 | CCL5 |
| CCND1 | CCND1 |

| | |
|----------------|--|
| CD27 | CD27 |
| CD274 | CD274 |
| CD276 | CD276 |
| CD3E | CD3E |
| CD4 | CD4 |
| CD40 | CD40 |
| CD40LG | CD40LG |
| CD44 | CD44 |
| CD47 | CD47 |
| CD68 | CD68 |
| CD74 | CD74 |
| CD86 | CD86 |
| CD8A | CD8A |
| CMKLR1 | CMKLR1 |
| CSF1R | CSF1R |
| CTLA4 | CTLA4 |
| CTNNB1 | CTNNB1 |
| CXCL10 | CXCL10 |
| CXCL9 | CXCL9 |
| CXCR6 | CXCR6 |
| DKK2 | DKK2 |
| EPCAM | EPCAM |
| FAS | FAS |
| FOXP3 | FOXP3 |
| GZMB | GZMB |
| HAVCR2 | HAVCR2 |
| HIF1A | HIF1A |
| HLA-DQ | HLA-DQA1 |
| HLA-DRB | HLA-DRB1,HLA-DRB3,HLA-DRB5,HLA-DRB4,HLA-DPB1 |
| HLA-E | HLA-E |
| ICAM1 | ICAM1 |
| ICOSLG | ICOSLG,LOC102723996 |
| IDO1 | IDO1 |
| IFNAR1 | IFNAR1 |
| IFNG | IFNG |
| IFNGR1 | IFNGR1 |
| IL12B | IL12B |
| IL15 | IL15 |
| IL6 | IL6 |
| ITGAM | ITGAM |
| ITGAV | ITGAV |
| ITGAX | ITGAX |
| ITGB2 | ITGB2 |
| ITGB8 | ITGB8 |

| | |
|---------------------------------|--|
| KRT | KRT18,KRT10,KRT17,KRT19,KRT6A,KRT6B,KRT6C,KRT7, KRT14 |
| LAG3 | LAG3 |
| LY6E | LY6E |
| MKI67 | MKI67 |
| MS4A1 | MS4A1 |
| NKG7 | NKG7 |
| pan-melanocyte | PMEL,S100B,SOX10 |
| PDCD1 | PDCD1 |
| PDCD1LG2 | PDCD1LG2 |
| PECAM1 | PECAM1 |
| PSMB10 | PSMB10 |
| PTEN | PTENP1,PTEN |
| PTPRC | PTPRC |
| STAT1 | STAT1 |
| STAT2 | STAT2 |
| STAT3 | STAT3 |
| TBX21 | TBX21 |
| TIGIT | TIGIT |
| TNF | TNF |
| TNFRSF9 | TNFRSF9 |
| VEGFA | VEGFA |
| VSIR | VSIR |
| Internal Reference Genes | |
| OAZ1 | OAZ1 |
| POLR2A | POLR2A |
| RAB7A | RAB7A |
| SDHA | SDHA,SDHAP1,SDHAP2,LOC220729,SDHAP3 |
| UBB | UBB,UBBP4 |
| Negative Controls | |
| NegPrb1 | |
| NegPrb2 | |
| NegPrb3 | |
| NegPrb4 | |
| NegPrb5 | |
| NegPrb6 | |

2.2 Cancer Transcriptomic Atlas (CTA)/Next Generation Sequencer: +1,800 genes

| Signaling pathways covered by CTA | | |
|--|--|--|
| Table 1. Adaptive Immunity Adaptive Immunity B cells BCR Signaling | Table 4. Cell Function Cell Function Apoptosis Autophagy | Table 7. Signaling Pathways Signaling Pathways AMPK Signaling Androgen Signaling |

| | | |
|-----------------------------------|--|-------------------------------------|
| Cancer Antigens | Cell Adhesion & Motility | EGFR Signaling |
| MHC Class I Antigen Presentation | Cell Cycle | ERBB2 Signaling |
| MHC Class II Antigen Presentation | Cilium Assembly | Estrogen Signaling |
| T-cell Checkpoints | Differentiation | FGFR Signaling |
| T cells | DNA Damage Repair | FoxO Signaling |
| TCR Signaling | EMT | GPCR Signaling |
| TH1 Differentiation | Endocytosis | Hedgehog Signaling |
| TH17 Differentiation | Epigenetic Modification | HIF1 Signaling |
| TH2 Differentiation | Immortality & Stemness | Insulin Signaling |
| TH9 Differentiation | Ion Transport | JAK-STAT Signaling |
| Treg Differentiation | Lysosome | MAPK Signaling |
| Table 2. Innate Immunity | Oxidative Stress | MET Signaling |
| Innate Immunity | Phagocytosis | mTOR Signaling |
| Complement System | Protectoxic Stress | Myc |
| Dendritic Cells | RNA Processing | NO Signaling |
| DNA Sensing | Senescence | Notch Signaling |
| Glycan Sensing | Table 5. Metabolism | p53 Signaling |
| Host Defense Peptides | Metabolism | PDGF Signaling |
| Inflammasomes | Amino Acid Synthesis & Transport | PI3K-Akt Signaling |
| Myeloid Inflammation | Arginine Metabolism | PPAR Signaling |
| Neutrophil Degranulation | Fatty Acid Oxidation | Purinergic Signaling |
| NK Activity | Fatty Acid Synthesis | Retinoic Acid Signaling |
| NLR Signaling | Glutamine Metabolism | TGF-beta Signaling |
| RAGE Signaling | Glycolysis & Glucose Transport | VEGF Signaling |
| RNA Sensing | Glycosylation | Wnt Signaling |
| TLR Signaling | IDH1/2 | Table 8. Tissue Compartment |
| Table 3. Immune Response | Lipid Metabolism | Tissue Compartment |
| Immune Response | Mitochondrial Metabolism / TCA | Tumor |
| Chemokine Signaling | Nucleotide Synthesis | Immune |
| Cytotoxicity | Pentose Phosphate | Stroma |
| IL-1 Signaling | Tryptophan & Kynurenine Metabolism | Table 9. Biological Category |
| IL-17 Signaling | Vitamin & Cofactor Metabolism | Tissue Compartment |
| IL-2 Signaling | Table 6. Physiology & Disease | Tumor Biology |
| IL-6 Signaling | Physiology & Disease | Immune Response |
| Immune Exhaustion | Angiotensin System | Microenvironment |
| Interferon Response Genes | Circadian Clock | |
| Lymphocyte Regulation | Drug Resistance | |
| Lymphocyte Trafficking | Glioma | |
| NF-kB Signaling | Leukemia | |
| Other Interleukin Signaling | Matrix Remodeling and Metastasis | |
| Prostaglandin Inflammation | Melanoma | |
| TNF Signaling | Neuroendocrine Function | |
| Type I Interferon Signaling | Prostate Cancer | |
| Type II Interferon Signaling | | |

Type III Interferon Signaling

-----END-----

Annex 3. Checklist for selection of cases and tissue qualification for DSP assay
GeoMX DSP check list for selection of cases and tissue qualification
DSP project #

1. Pathology assessment of samples:

Check if the following steps have been completed:

H&E slides scanned in Aperio System.

-Aperio Project ID: _____

File with Clinicopathological information received

Check if the following information was received:

MRN

Tissue ID

Time Point. If the patient has two TPs (TP0, TP1), check the experimental design for specific instructions on sample selection.

Final Diagnosis. (Apollo provides with Primary and Reviewed diagnosis).

Organ site.

Type of biopsy (surgical specimen, core biopsy, endoscopic biopsy, TMA) (A map of the TMA must be requested)

Localization (primary organ vs metastasis in lymph nodes or other organs).

Other

2. Quality control:

Check if the following information was evaluated

Presence of tumor cells (Yes/No): > than 100 viable tumor cells.

Presence of normal adjacent tissue.

Presence of Tumor associated Tertiary lymphoid aggregates (TLS).

Presence of Invasive margin (IM).

% of tumor infiltrating lymphocytes (TILs)*.

Comments: Annotate potential presence of artifacts that can interfere with the staining (background) including blood cells, fibrin, bone or necrosis. Areas with artifacts, mucin or necrosis as well as samples with inadequate processing issues will be excluded for the further selection of regions of interest (ROIs).

3. Tissue qualification for DSP assay:

Check if the following steps have been completed:

Sample selected for DSP assay annotated.

ROI plan annotated : estimated number of ROI according to research strategy and type of biopsy.

AOI plan annotated: estimated number of AOI based in segmentation strategy and ROI.

4. Requisition of material to perform DSP assay:

Check if the following steps have been completed:

Unstained sections requested (Date)

5. Storage of slides:

Check if the following steps have been completed:

Unstained slides received.

Date:

Location:

Total number of sections: _____

-Number of sections per sample/case

Total number of sample/cases:

6. Quality assurance of unstained sections:

Check if the following steps have been completed:

Unstained sections are similar to H&E previously reviewed?

-Yes.

-No.

Note: If unstained section varies significantly to the H&E slide reviewed, an H&E of a serial section should be requested for optimal ROI selection.

Number of slides with adequate tissue sections: _____

(Tissue sections do not have excessive folding or detachment. Tissue section fit the DSP scanning area. There is not tissue in the tip calibration area)

Slide with issues: _____

__ Tissue sections are suboptimal. DSP can be performed with limitations:

Slide ID:

__ Tissue sections are not adequate. DSP cannot be performed:

Slide ID:

The following issues have being observed:

__ Tissue sections have excessive folding

__ Tissue sections have extensive detachment.

__ Tissue sections do not fit the DSP scanning area.

__ There is tissue in the tip calibration area.

__ Other

__ unstained section selected for DSP assay.


-----END-----

Annex 4: Selection of ROI in the scanned images using GeoMx Digital image analysis software

Selection of ROI and segmentation using the GeoMx Digital image analysis software

1. Prepare pre-assessment files including excel file containing information of samples and H&E images with the potential placement of ROI. H&E image slides can be visualized using Aperio Image Scope (eSlide manager) or Halo (IndicaLab). If a grid with Tumor Infiltrating Lymphocytes (TILs) % was prepared, open file containing the information.
2. Log in into GeoMx DSP software through Google Chrome based app
<https://10.104.100.32/login> (DSP#1)
<https://10.104.100.33/login> (DSP#2)
3. Search for the scanned slide at the bottom of the screen (from slot 1 to slot 4), and click on the selected slide. The image should consist of immune fluorescent slide with multiple channels.
4. Review that the DSP ID (located at the top left of the image) and scan number information (S1;, S2, etc.) in the image slide match the corresponding ID from file.
5. Adjust intensity of the IF staining by scrolling to desire channel and adjusting threshold. Pick desired color for the selected channer image. Compare the scanned tissue with the H&E slide for orientation (screenshot 1). The laboratory should have provided information on the morphology markers, fluorophores and DSP channel.

Screenshot 1



6. A quick guide for DSP channels and fluorophores is provided below (screenshot 2).
 Screenshot 2:

| DSP channel / Em max | Excitation range | Emission range | Emission range |
|----------------------|------------------|----------------|---|
| FITC / 525 nm | 480 +/- 28 nm | 516 +/- 23 nm | AF488*, FITC, DyLight488, Fluorescein, 5-FAM, Rhodamine Green |
| Cy3 / 568 nm | 538 +/-19 nm | 564 +/-15 nm | AF532*, Cy3, PE, TRITC |
| Texas Red / 615 nm | 588 +/- 19 nm | 623 +/- 30 nm | AF594*, Texas Red, DyLight594 |
| Cy5 / 666 nm | 645 +/- 19 nm | 683 +/- 30 nm | AF647*, DyLight650, Cy5, AF660w |

* recommended first choice




FIGURE 2: Fluorescence SpectraViewer (ThermoFisher) displaying DSP's excitation and emission filters as well as excitation (dotted lines) and emission (solid, filled in graphs) for the four recommended fluorophores. AF488 = green, AF532 = yellow, AF594 = orange, AF647 = red.

7. The information of any specific slide can be accessed by clicking in “scan parameters” (screenshot 3)

Screenshot 3:

Channel Settings

| | Channel / Emission Max | Fluorophore | Biological Target | Biological Class | Exposure Time |
|-------------------------------------|------------------------|-------------|-------------------|------------------|---------------|
| <input checked="" type="checkbox"/> | FITC/525 nm | SYTO 13 | DNA | | 50 ms |
| <input checked="" type="checkbox"/> | Cy3/568 nm | Alexa 532 | PanCK | | 400 ms |
| <input checked="" type="checkbox"/> | Texas Red/615 nm | Alexa 594 | CD45 | | 400 ms |
| <input checked="" type="checkbox"/> | Cy5/666 nm | Cy5 | CD 68 | | 400 ms |

8. Select all ROIs in the sample following the experimental design strategy. Note that areas with artifacts should be excluded.
9. Generate segments (e.g. Tumor; Stroma; and Immune) according to the experimental design strategy and replicate the same strategy for segment generation including order and segment definitions in each sample from the same experiment to keep consistency.
 - Depending of the type of sample, additional segments can be generated to exclude background or empty spaces based on biomarker expression patterns (e.g. Glass; Elastic fibers; Blood).
 - Pick a desired color per segment; we recommend using different colors to segment mark-ups to have a better contrast with the visualization markers.
10. Click in advanced parameters for segment generation: Choose desired Hole size and Particle size settings (The initial settings recommended for tumor segment is Hole size:1000 and Particle size:50; and for Stroma or Immune segments Hole size: 160 Particle size: 10 (immune) or 50 (stroma); Erode um: 1.). Settings can be tailored to specific projects based in biomarker and architectural patterns of biological compartments).
11. Save the changes.
12. Adjust immunofluorescence intensities thresholds by each ROI to reach the best segmentation parameters. Turn on and off the mark-ups to make sure segmentation parameters are optimal.
13. Note that changing parameters for segmentation after this step will reset IF intensities and thresholds should be adjusted accordingly.
14. Example of Segment definitions and channel thresholds is shown in screenshot 4:

Screenshot 4:

The screenshot displays the 'Segments' panel in the nanoString GeoMx DSP software. It includes a 'Segment Definitions' section with two defined segments: 'Tumor' (cyan) and 'Stroma' (yellow). The 'Tumor' segment is defined with parameters: Erode: 1 μm, N-Dilate: 2 μm, Hole Size: 1000 μm², Particle Size: 50 μm², and Collection Order: 1. The 'Stroma' segment is defined with parameters: Erode: 1 μm, N-Dilate: 2 μm, Hole Size: 160 μm², Particle Size: 50 μm², and Collection Order: 2. Below this is a 'Channel Thresholds' table.

| RO1 | FITC 525 nm | Cy3 568 nm | Texas Red 615 nm | Cy5 666 nm |
|--------------------------|----------------|---------------|---------------------|---------------|
| 001 | 5 | 18 | n/a | n/a |
| 002 | 5 | 18 | n/a | n/a |
| 003 | 5 | 21 | n/a | n/a |
| 004 | 5 | 13 | n/a | n/a |
| 005 | 5 | 24 | n/a | n/a |
| 006 | 5 | 10 | n/a | n/a |
| 007 | 1 | 18 | n/a | n/a |
| 008 | 5 | 10 | n/a | n/a |
| 009 | 5 | 23 | n/a | n/a |
| 010 | 5 | 12 | n/a | n/a |
| 011 | 5 | 12 | n/a | n/a |
| 012 | 1 | 20 | n/a | n/a |
| 013 | 1 | 13 | n/a | n/a |
| Override All (if > 0) | 0 | 0 | n/a | n/a |

The main image shows a histological section with the 'Tumor' and 'Stroma' segments highlighted in cyan and yellow, respectively. The background is black, and there are other colored regions in blue and red.

15. Exclude any segment generated that do not qualify for DSP analysis including segments with suboptimal number of cells and segments generated to avoid artifacts (e.g. glass, elastic fibers)
16. Review the count of nuclei by segment and exclude any segment with less than 20 nuclei for protein and less than 200 cells for RNA. DSP counts of nuclei may show an inaccurate nuclei count, then it is recommended to confirm number of nuclei by direct observation of image (screenshot 5).

Screenshot 5

The screenshot displays the 'ROIs' section of the software. At the top, there is a blue header with 'ROIs' and a back arrow. Below this, a message states: 'This scan is not currently editable. You can see any Regions Of Interest (ROIs) that were selected for it below.' There are two buttons: 'Export All' and 'Hide Segments'. Below these are 'Export' buttons for each ROI.

ROI: 006
 Rectangle - Center: (**5347.0** , **3574.8**) Size: (**660.0** x **785.0** μm)
 Area: 518100.0 μm^2 , Nuclei Count: 3221
 Created: 05/17/2022, 12:49:29 PM

| Segments: | Name | Area (μm^2) | Nuclei | |
|-----------|-----------|--------------------------|--------|--|
| | 1) Tumor | 274273.5 | 2206 | |
| | 2) Stroma | 213533.8 | 1069 | |

ROI: 007
 Rectangle - Center: (**5380.9** , **1874.3**) Size: (**660.0** x **785.0** μm)
 Area: 518100.0 μm^2 , Nuclei Count: 2961
 Created: 05/17/2022, 12:49:36 PM

| Segments: | Name | Area (μm^2) | Nuclei | |
|-----------|-----------|--------------------------|--------|--|
| | 1) Tumor | 154892.3 | 1363 | |
| | 2) Stroma | 313264.9 | 1636 | |

17. Once this is finalized and segments have passed pathology QC, send AOI for collection.
18. Continue the selection of ROI and segmentation in the next slide.
19. Communicate to the lab when Selection of ROI have been completed, so the lab can proceed with hybridization and proper steps for mapping back the DSP counts to image slide.
20. After the DSP counts are mapped back to DSP device. An excel file should be ready in order to merge pathological annotation details to the GeoMx DSP analysis suite software. This file can be filled in while performing segmentation or after the DSP counts are mapped back in order to confirm ROI and segment information.
21. The excel file should contain the following information.
 - a. Batch_HE: *Batch number of H&E slides received*
 - b. Slide_ID_DSP: *This Slide ID must match the Slide ID used to label the slides scanned in the DSP device. Of note, a second scan or repeated cases of the same sample need to be coded with "t" or "s".*
 - i. "t" if new sections from the same block is stained and scanned.
 - ii. "s" if new scan from the same stained slide is performed (e.g. PA03, PA03s2)
 - c. ROI (label): *This cell must match ROI (label) in the DSP annotated regions. It needs to be entered with apostrophe (e.g. '001) to match DSP template file*
 - d. Segment (Name / Label): *This cell must match Segment label used when segmenting images in the DSP device.*
 - e. Case_ID: *Unique identification of the sample included in the analysis.*
 - f. TID: *Tissue ID provided by collaborators.*

- g. Primary_Dx: *Information extracted from collaborators.*
- h. Strategy: *Indicates the strategy used in each ROI (e.g. Tumor vs. Stroma; tumor vs. Immune, no segmentation).*
- i. ROI_Type: *e.g. Intratumoral vs Invasive margin vs. TLS per experimental design strategy.*
- j. Compartment: *Biological compartment after segmentation or the whole ROI if not segmented (e.g. tumor vs. Immune).*
- k. QC_CK: *Quality of staining of tumor biomarker (e.g. Cytokeratin) (Ok, background; low).*
- l. QC_CD45: *Quality of staining of CD45 in immune cells (Ok, background; low).*
- m. QC_Number_nuclei_Immune: *Quality control of the number of nuclei obtained from the "Immune" segment by DSP platform (visual assesment; e.g. ok/hypersegmented/hiposegmented).*
- n. QC_Number_nuclei_Tumor: *Quality control of the number of nuclei obtained from the "Tumor" segment by DSP platform (visual assesment; e.g. ok/hypersegmented/hiposegmented).*
- o. QC_Number_nuclei_Stroma: *Quality control of the number of nuclei obtained from the "Stroma" segment by DSP platform (visual assesment; e.g. ok/hypersegmented/hiposegmented).*

-----END-----