



16S rRNA with MiSeq

Author: Gerold Bongers
SOP: MTC004
Revision: 1
Effective date: 30 Nov 2018

1 Abstract

The Purpose of this protocol is to define the steps for the preparation and sequencing of 16S rRNA gene sequence libraries using the Illumina MiSeq sequencing platform. Method adapted from Kozich, J. J. et al Appl. Environ. Microbiol. 79, 5112–20 (2013) [https://github.com/SchlossLab/MiSeq_WetLab_SOP/blob/master/MiSeq_WetLab_SOP_v4.md]

2 Materials and equipment

2.1 Reagents

- Hard-Shell® 96-Well PCR Plates, low profile, thin wall, skirted, white/clear [Bio-Rad, HSP9601]
- PCR Plate Heat Seal, foil, peelable [Bio-Rad, 1814045]
- UltraPure™ DNase/RNase-Free Distilled Water [Life Technologies, 10977015]
- Phusion High-Fidelity PCR Master Mix [NEB, M0531L]
- Quant-iT™ dsDNA Assay Kit, broad range [Life Technologies, Q33130]
- Corning™ 96-Well Solid Black Plates [Fisher Scientific, 07-200-590]
- Beckman Coulter AMPURE XP [Fisher Scientific, NC9959336]
- Misc (tips, gel)
- Biomek tips [Beckman Coulter, Biomek AP96 P250 & Span P250]
- Primers [IDT, NA]
- ZymoBIOMICS Microbial Community Standard [Zymo Research, D6300]
- TE, pH 8.0 [Life Technologies, AM9849]
- Buffer EB [Qiagen, 19086]
- Quarter Reservoir Sterile Divided by Length [Beckman Coulter, 372788]
- GEL LOAD DYE PURP (6X) 4 ML [Fisher Scientific, 50591186]
- QUICK-L 2-LOG DNA LADDER [Fisher Scientific, 50994896]
- SYBR SAFE DNA GEL STAIN [Fisher Scientific, S33102]
- AGAROSE LOW EEO 100G [Fisher Scientific, BP160100]
- NaCl (5 M) [Life Technologies, AM9760G]
- PhiX Control v3 [Illumina, C-110-3001]
- MiSeq Reagent Kit v2 (300-cycles) [Illumina, MS-102-2002]

2.2 Equipment

- Plate shaker [Benchmark, Orbishaker]

- Plate sealer [Bio-Rad, PCR Plate Sealer]
- PCR machine [Bio-Rad, C1000 Touch Thermal Cycler]
- Gel Doc [Bio-Rad, Gel Doc EZ Imager]
- Microplate reader, filter 450 nm [BioTek, Synergy HTX]
- Liquid Handler [Beckman Coulter, Biomek FXP]
- Microcentrifuge (>3000g) [Eppendorf, 5424]
- Gel cast [Thermo Scientific, Easycast B1A, B2]
- Fluorometer [ThermoFisher, Qubit® 2.0 Fluorometer]

3 Safety precautions

3.1 Human specimen handling

Universal precautions should be followed when handling human samples considered potentially infectious. Minimal protective clothing requirements are lab coat, gloves and safety glasses.

3.2 Waste

Dispose of all waste as biohazardous material.

3.3 Specific precautions

4 Method

4.1 Considerations

4.1.1 Low biomass samples

To increase the amount of input DNA for low microbial biomass samples (run a test PCR to optimize against a no-template control). Consider using a using DNA isolation kits designed to minimize kit contamination (e.g., the QIAamp UCP (UltraClean production) Pathogen Mini Kit (QIAGEN)). Consider pooling samples and run of the same extraction column to increase microbial biomass. Work in a down-flow cabinet to prevent environmental contamination.

4.1.2 Controls

DNA extraction controls: on each plate include ZymoBIOMICS Microbial Community Standard, 3 longitudinal fecal reference samples obtained from healthy volunteers and a no sample control on random positions. PCR controls: no template control and Clemente lab gut reference sample.

4.2 Setup

4.2.1 Primer dilution

1. Reconstitute indexed primer stock to 100 μ M in TE (see Appendix for primer design).
2. Reconstitute sequence primers to 100 μ M in EB (see Appendix for primer design).
3. Dilute Sx5xx primers to 10 μ M (140 μ l stock + 1360 μ l EB), mix by vortexing.
4. Dilute Sx7xx primers to 10 μ M (100 μ l stock + 900 μ l EB), mix by vortexing.
5. Add 2x 720 μ l late for Sx5xx and 1x1000 μ l for Sx7xx to 96 deep well plate
6. With a multichannel transfer 50 μ l of Sx5xx and 50 μ l of Sx7xx to 4 PCR 96 well plates (See Appendix for plate layout) according to the scheme: Plate 1: A701 – A712 with A501 – A508; Plate 2: A701 – A712 with B501 – B508; Plate 3: B701 – B712 with B501 – B508; Plate 4: B701 – B712 with A501 – A508.


7. Seal plates for 3 sec @ 180C.
8. Vortex plates for 3 min at 300 rpm and spin at 1000g for 1 min.
9. Barcode plates and store at -80C

4.2.2 Dilute DNA

10. Quantify the DNA in each sample using the Quant-iT dsDNA Assay Kit, Broad Range (MTCP007-Quant-iT_dsDNA_Assay).

Note: If the concentration is less than 50ng/μL by the Qubit® dsDNA BR Assay Kit, the samples should be reanalyzed using the Qubit® dsDNA HS Assay Kit.

11. Standardize all samples to 2 ng/ul using BioMek [protocol: MTC/2_DiluteDNA]. Specify dilutions in template DNA_Dil.csv. Make sure to load the correct DNA_Dil.csv file for both steps and setup deck as per figure below. Fill the first well of the Modular Reservoir with EB. Home the axis before use under the instrument menu. Make sure all PCR plates are labeled/barcoded and take a picture of the deck before starting the protocol.



12. After the protocol completes, vortex plates for 3 min at 300 rpm and spin at 1000g for 1 min.
13. Seal plates for 3 sec @ 180C.
14. Store plates at -80C.

4.3 PCR

4.3.1 Master mix (96 samples)

15. Combined the following:

Master Mix: 96 samples		
Reagent	ul	Total
H2O	8.0	845
2xPhusion	12.5	1320


4.3.2 Reaction

16. Barcode PCR plates.
17. Transfer Master Mix with repeat pipette to PCR plate.
18. Transfer DNA with multichannel to PCR plate.
19. Transfer Paired Primers with multichannel to PCR plate.

Reagent	Unit
DNA (2ng/ul)	2.0 ul
Paired primers	2.5 ul
Master Mix	20.5 ul
Total	25 ul

20. Repeat for up to four 96 well plates. Seal plates, vortex briefly and spin down

Note: Alternatively use Biomek Protocol MTC/3_SetupPCR. Setup the deck as indicated below and take a picture of the deck before starting the protocol.



4.3.3 PCR program

21. Place in thermocycler and run the program below


Step	Temp	Time
Denaturation	98	0:30
25 cycles	98	0:10
	50	0:30
	72	0:30
Extension	72	2:00
Hold	4	Inf

4.4 Gel Electrophoresis

22. 1 or 2 random rows of 12 should be selected from each PCR plate and run on a gel to confirm success of the PCR.
23. Use 5 ul of sample, 1 ul of 6X loading dye in a 1% agarose gel (bring to boil 0.5g agarose in 50 ml 1xTAE, cool down and add 5 ul of 10.000x SYBR safe).
24. Run at 100v for 30 minutes alongside a 2-log ladder.
25. Photograph gel using the EZ Imager. Check to be sure there is a band for every well


4.5 Normalization and Pooling

26. Quantify the DNA in each sample using the Quant-IT dsDNA Assay Kit, Broad Range (MTCP007-Quant-iT_dsDNA_Assay).
27. Combine in ~equal proportions with Biomek (Protocol MTC/4_CombineLibrary) following the Library.csv template. Setup deck as below and take a picture before starting. Place a barcoded library Eppendorf tube in position A1 of the Small TubeRack. Transfer a maximum of 10ul for any sample including NTC controls.



28.

4.6 Double size selection and cleanup using AMPure XP Beads



4.6.1 Right side selection

This step is used to bind the large, unwanted fragments to the beads. The supernatant will contain the desired fragments.

29. Add 77 μ l (0.77X) resuspended AMPure XP beads to 100 μ l DNA solution. Mix well on a vortex mixer or by pipetting up and down at least 10 times.
30. Incubate for 5 minutes at room temperature.
31. Place the tube on a magnetic rack to separate the beads from the supernatant. After the solution is clear (about 5 minutes), carefully transfer 100 μ l of the supernatant to a new. Discard beads that contain the large fragments.

4.6.2 Left side selection

This step will bind the desired fragment sizes to the beads. Unwanted smaller fragment sizes will not bind to the beads.

32. Add 33 μ l (0.33X) resuspended AMPure XP beads to 100 μ l of the supernatant, mix well and incubate for 5 minutes at room temperature.
33. Put the tube on a magnetic rack to separate beads from supernatant. After the solution is clear (approximately 3 minutes), carefully remove and discard the supernatant. Be careful not to disturb the beads that contain DNA targets (Caution: do not discard beads).
34. Add 200 μ l of 80% freshly prepared ethanol to the tube while in the magnetic stand. Incubate at room temperature for 30 seconds, and then carefully remove and discard the supernatant.
35. Repeat Step above once.
36. Keeping the tube on the magnetic rack, with the cap open, air dry the beads for 5 minutes.
Caution: Do not overdry the beads. This may result in lower recovery of DNA target.
37. Remove the tube from the magnet. Elute DNA target from beads into 42 μ l sterile 0.1X TE. Mix well on a vortex mixer or by pipetting up and down, incubate for 2 minutes at room temperature.
38. Put the tube in a magnetic rack until the solution is clear, approximately 3 minutes.
39. Transfer approximately 40 μ l of the supernatant to a barcoded clean tube.

4.7 Gel Electrophoresis

40. Use 5 μ l of sample (before and after cleanup), 1 μ l of 6X loading dye in a 1% agarose gel (bring to boil 0.5g agarose in 50 ml 1xTAE, cool down and add 5 μ l of 10.000x SYBR safe).
41. Run at 100v for 30 minutes alongside a 2-log ladder.
42. Photograph gel using the EZ Imager. Compare before and after AMPure purification ensure that there is a single band after cleanup. Alternatively, use a BioAnalyzer or DNA ScreenTape to evaluate Library.

4.8 Dilution

43. Quantify Library sample (Qubit HS, follow manufactures protocol) and normalize to 15 nM using the following formula (library size V4 386 bp).

$$\frac{(\text{concentration in ng}/\mu\text{l})}{(660 \text{ g/mol} \times \text{average library size in bp})} \times 10^6 = \text{concentration in nM}$$

44. Barcode tube and store library at -80C or send for MiSeq 2x250bp sequencing.

5 Appendix

5.1 Primers

5.1.1 Overall design considerations:

- The sequencing primers must have a melting temperature near 65°C. This can be achieved by altering the pad sequence.
- The index sequences must balance the number of bases at each position. The index sequences listed here have a 25% ATGC composition at each site. If you are going to cherry pick indices from the list, make sure that you have even representation.

5.1.2 Generic PCR primer design:

AATGATACGGCGACCACCGAGATCTACAC <i5><pad><link><16Sf> VX.N5??

CAAGCAGAAGACGGCATACGAGAT <i7><pad><link><16Sr> VX.N7??

Generic read 1 primer design: <pad><link><16Sf> VX.read1

Generic read 2 primer design: <pad><link><16Sr> VX.read2

Generic index read primer design: Reverse complement of (<pad><link><16Sr>) VX.p7_index

The listed sequences in the generic design, above, are the adapter sequences to allow annealing of the amplicons to the flow cell. The i5 and i7 sequences are the 8-nt index sequences. The pad is a 10-nt sequence to boost the sequencing primer melting temperatures. The link is a 2-nt sequence that is anti-complementary to the known sequences. The 16Sf and 16Sr are the gene specific primer sequences. Primers are purchased from IDT with no special purification. This system should work for any other region of the 16S rRNA gene or any other gene. The only thing to change would be the 16Sf/16Sr sequences and confirm that when combined with the pad sequence that the melting temperature is near 65°C.

5.1.3 Sequence primers

Name	pad	link	16Sf
Read1	TATGGTAATT	GT	GTGCCAGCMGCCGCGTAA
Name	pad	link	16Sr
Read2	AGTCAGTCAG	CC	GGACTACHVGGGTWCTAAT
Name	pad	link	16Sr
Readl	AGTCAGTCAG	CC	GGACTACHVGGGTWCTAAT

5.1.4 Forward index primers

Name	adapter	i5	pad	link	16Sf
SA501	AATGATACGGCGACCACCGAGATCTACAC	ATCGTACG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SA502	AATGATACGGCGACCACCGAGATCTACAC	ACTATCTG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SA503	AATGATACGGCGACCACCGAGATCTACAC	TAGCGAGT	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SA504	AATGATACGGCGACCACCGAGATCTACAC	CTGCGTGT	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SA505	AATGATACGGCGACCACCGAGATCTACAC	TCATCGAG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SA506	AATGATACGGCGACCACCGAGATCTACAC	CGTGAGTG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SA507	AATGATACGGCGACCACCGAGATCTACAC	GGATATCT	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SA508	AATGATACGGCGACCACCGAGATCTACAC	GACACCGT	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SB501	AATGATACGGCGACCACCGAGATCTACAC	CTACTATA	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SB502	AATGATACGGCGACCACCGAGATCTACAC	CGTTACTA	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SB503	AATGATACGGCGACCACCGAGATCTACAC	AGAGTCAC	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SB504	AATGATACGGCGACCACCGAGATCTACAC	TACGAGAC	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SB505	AATGATACGGCGACCACCGAGATCTACAC	ACGTCTCG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SB506	AATGATACGGCGACCACCGAGATCTACAC	TCGACGAG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SB507	AATGATACGGCGACCACCGAGATCTACAC	GATCGTGT	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SB508	AATGATACGGCGACCACCGAGATCTACAC	GTCAGATA	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SC501	AATGATACGGCGACCACCGAGATCTACAC	ACGACGTG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SC502	AATGATACGGCGACCACCGAGATCTACAC	ATATACAC	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SC503	AATGATACGGCGACCACCGAGATCTACAC	CGTCGCTA	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SC504	AATGATACGGCGACCACCGAGATCTACAC	CTAGAGCT	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SC505	AATGATACGGCGACCACCGAGATCTACAC	GCTCTAGT	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SC506	AATGATACGGCGACCACCGAGATCTACAC	GACACTGA	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SC507	AATGATACGGCGACCACCGAGATCTACAC	TGCGTACG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SC508	AATGATACGGCGACCACCGAGATCTACAC	TAGTGTAG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SD501	AATGATACGGCGACCACCGAGATCTACAC	AAGCAGCA	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SD502	AATGATACGGCGACCACCGAGATCTACAC	ACGCGTGA	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SD503	AATGATACGGCGACCACCGAGATCTACAC	CGATCTAC	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SD504	AATGATACGGCGACCACCGAGATCTACAC	TGCGTCAC	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SD505	AATGATACGGCGACCACCGAGATCTACAC	GTCTAGTG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SD506	AATGATACGGCGACCACCGAGATCTACAC	CTAGTATG	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SD507	AATGATACGGCGACCACCGAGATCTACAC	GATAGCGT	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA
SD508	AATGATACGGCGACCACCGAGATCTACAC	TCTACT	TATGGTAATT	GT	GTGCCAGCMGCCGCGGTAA

5.1.5 Reverse index primers

Name	adapter	i7	pad	link	16Sr
SA701	CAAGCAGAAGACGGCATAACGAGAT	AACTCTCG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA702	CAAGCAGAAGACGGCATAACGAGAT	ACTATGTC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA703	CAAGCAGAAGACGGCATAACGAGAT	AGTAGCGT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA704	CAAGCAGAAGACGGCATAACGAGAT	CAGTGAGT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA705	CAAGCAGAAGACGGCATAACGAGAT	CGTACTCA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA706	CAAGCAGAAGACGGCATAACGAGAT	CTACGCAG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA707	CAAGCAGAAGACGGCATAACGAGAT	GGAGACTA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA708	CAAGCAGAAGACGGCATAACGAGAT	GTCGCTCG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA709	CAAGCAGAAGACGGCATAACGAGAT	GTCGTAGT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA710	CAAGCAGAAGACGGCATAACGAGAT	TAGCAGAC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA711	CAAGCAGAAGACGGCATAACGAGAT	TCATAGAC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SA712	CAAGCAGAAGACGGCATAACGAGAT	TCGCTATA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB701	CAAGCAGAAGACGGCATAACGAGAT	AAGTCGAG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB702	CAAGCAGAAGACGGCATAACGAGAT	ATACTTCG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB703	CAAGCAGAAGACGGCATAACGAGAT	AGCTGCTA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB704	CAAGCAGAAGACGGCATAACGAGAT	CATAGAGA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB705	CAAGCAGAAGACGGCATAACGAGAT	CGTAGATC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB706	CAAGCAGAAGACGGCATAACGAGAT	CTCGTTAC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB707	CAAGCAGAAGACGGCATAACGAGAT	GCGCACGT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB708	CAAGCAGAAGACGGCATAACGAGAT	GGTACTAT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB709	CAAGCAGAAGACGGCATAACGAGAT	GTATACGC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB710	CAAGCAGAAGACGGCATAACGAGAT	TACGAGCA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB711	CAAGCAGAAGACGGCATAACGAGAT	TCAGCGTT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SB712	CAAGCAGAAGACGGCATAACGAGAT	TCGCTACG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC701	CAAGCAGAAGACGGCATAACGAGAT	ACCTACTG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC702	CAAGCAGAAGACGGCATAACGAGAT	AGCGCTAT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC703	CAAGCAGAAGACGGCATAACGAGAT	AGTCTAGA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC704	CAAGCAGAAGACGGCATAACGAGAT	CATGAGGA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC705	CAAGCAGAAGACGGCATAACGAGAT	CTAGCTCG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC706	CAAGCAGAAGACGGCATAACGAGAT	CTCTAGAG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC707	CAAGCAGAAGACGGCATAACGAGAT	GAGCTCAT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC708	CAAGCAGAAGACGGCATAACGAGAT	GGTATGCT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC709	CAAGCAGAAGACGGCATAACGAGAT	GTATGACG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC710	CAAGCAGAAGACGGCATAACGAGAT	TAGACTGA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC711	CAAGCAGAAGACGGCATAACGAGAT	TCACGATG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SC712	CAAGCAGAAGACGGCATAACGAGAT	TCGAGCTC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD701	CAAGCAGAAGACGGCATAACGAGAT	ACCTAGTA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD702	CAAGCAGAAGACGGCATAACGAGAT	ACGTACGT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD703	CAAGCAGAAGACGGCATAACGAGAT	ATATCGCG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD704	CAAGCAGAAGACGGCATAACGAGAT	CACGATAG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD705	CAAGCAGAAGACGGCATAACGAGAT	CGTATCGC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD706	CAAGCAGAAGACGGCATAACGAGAT	CTGCGACT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD707	CAAGCAGAAGACGGCATAACGAGAT	GCTGTAAC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD708	CAAGCAGAAGACGGCATAACGAGAT	GGACGTTA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD709	CAAGCAGAAGACGGCATAACGAGAT	GGTCGTAG	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD710	CAAGCAGAAGACGGCATAACGAGAT	TAAGTCTC	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD711	CAAGCAGAAGACGGCATAACGAGAT	TACACAGT	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT
SD712	CAAGCAGAAGACGGCATAACGAGAT	TTGACGCA	AGTCAGTCAG	CC	GGACTACHVGGGTWTCTAAT

5.1.6 Primer plate layout

		SA701	SA702	SA703	SA704	SA705	SA706	SA707	SA708	SA709	SA710	SA711	SA712
	Plate 1	1	2	3	4	5	6	7	8	9	10	11	12
SA501	A												
SA502	B												
SA503	C												
SA504	D												
SA505	E												
SA506	F												
SA507	G												
SA508	H												
		SA701	SA702	SA703	SA704	SA705	SA706	SA707	SA708	SA709	SA710	SA711	SA712
	Plate 2	1	2	3	4	5	6	7	8	9	10	11	12
SB501	A												
SB502	B												
SB503	C												
SB504	D												
SB505	E												
SB506	F												
SB507	G												
SB508	H												
		SB701	SB702	SB703	SB704	SB705	SB706	SB707	SB708	SB709	SB710	SB711	SB712
	Plate 3	1	2	3	4	5	6	7	8	9	10	11	12
SA501	A												
SA502	B												
SA503	C												
SA504	D												
SA505	E												
SA506	F												
SA507	G												
SA508	H												
		SB701	SB702	SB703	SB704	SB705	SB706	SB707	SB708	SB709	SB710	SB711	SB712
	Plate 4	1	2	3	4	5	6	7	8	9	10	11	12
SB501	A												
SB502	B												
SB503	C												
SB504	D												
SB505	E												
SB506	F												
SB507	G												
SB508	H												