

## SOP- 'V'/2017 - AQIX® RS-I Biospecimen Sample Procurement

**Preamble:** The following directives serve to indicate the versatility of the AQIX® RS-I solution in facilitating the procurement of samples of 'normal' or 'malignant' biospecimen tissue under either hypothermic or ambient storage/transport temperatures and the option of dividing the same sample into three segments allowing the operator to simultaneously process 'fresh', formalin fixed (FF) and snap frozen (SF) tissue samples over 1 - 72 hour procurement periods. Likewise, liquid biopsied cell samples (e.g., CTC's) may be processed in a similar manner.

### Method 1 – Hypothermic Storage/Transportation

- a) Check that the 8 month shelf-life period on the AQIX® RS-I container label has not expired.
- b) Store 125 mL containers of AQIX® RS-I solution @ **3 – 8 °C** in dark conditions.
- c) Transport AQIX® RS-I container @ 0-4 °C (over 'wet' ice) to sample retrieval site.
- d) Maintain AQIX® RS-I container @ 0-4 °C (over 'wet' ice) at sample retrieval site under dark conditions.
- e) Open a 125mL AQIX® RS-I container for the minimal time possible before inserting the sample.
- f) Quickly attach the lid closure and seal the 125 mL AQIX® RS-I container tightly.
- g) Transport back to laboratory @ 0-4 °C (over 'wet' ice) in a polystyrene, outer box.
- h) Process as a 'fresh' tissue sample using current proteomic, genomic, metabolomic and/or histochemical technologies, **or**,
- i) Further divide the sample into three segments, retain one 'fresh', snap freeze (SF) the second segment and formalin fix (FFPE) the last segment – process each segment accordingly, **or**,
- j) Store the AQIX® RS-I container + sample @ **3 – 8 °C** for up to 72 hours under dark conditions prior to proceeding as in h) or i) above.

### Method 2 – Ambient Storage/Transportation

- a) Check that the 8 month shelf-life period on the AQIX® RS-I container label has not expired
- b) Store 125 mL kit containers of AQIX® RS-I solution @ **15 - 25 °C** under dark conditions.
- c) Transport AQIX® RS-I container @ < 25 °C to sample retrieval site.
- d) Maintain AQIX® RS-I container @ < 25 °C under dark conditions at sample retrieval site.
- e) Open a 125mL AQIX® RS-I container for the minimal time possible before inserting the sample.
- f) Quickly attach the lid closure and seal the 125 mL container tightly.
- g) Transport back to laboratory @ < 25 °C in a polystyrene, outer box.
- h) Process as a 'fresh' tissue sample using current proteomic, genomic, metabolomic and/or histochemical technologies, **or**,
- i) Further divide the sample into three segments, retain one 'fresh', snap freeze (SF) the second segment and formalin fix (FFPE) the last segment – process each segment accordingly, **or**,
- j) Store the AQIX® RS-I container + sample @ **3 – 8 °C** for up to 72 hours under dark conditions prior to proceeding as in h) or i) above.

## Addendum

1. If tissue biopsy contamination becomes a problem then additional AQIX® RS-I solution may be needed to thoroughly rinse the biopsy sample before inserting into the 125 mL AQIX® RS-I container. Additionally, 50-100 mg/L of Chloromycetin or 25-50 mg/L Nanomycopolitin may be added to AQIX® RS-I solution to prevent bacterial contamination without compromising the viability of the biopsy sample.
2. The Table 'V' below indicates the versatility of the AQIX® RS-I solution in allowing the same liquid or solid biopsy sample to be subdivided into three and simultaneously preserved as a 'fresh', frozen (SF) and formalin fixed (FFPE) samples for diagnostic analysis.

<b>Table 'V' - AQIX® Technology – Versatile Tissue Biopsy Procurement Methods</b>		
<b>Tissue Collection to Prognostics</b>	<b>Current Practices</b>	<b>AQIX® RS-I Protocols</b>
<b>Collection and immediate holding</b>	Biopsy delivered dry/soaked in saline/ wrapped in gauze /formalin fixed/	<b>Biopsies simply immersed in requisite volume of RS-I @ 3 - 8°C and held or transported before initial pathologist examination</b>
<b>Transportation</b>	Snap frozen ( <b>SF</b> ), Formalin fixed ( <b>FFPE</b> ) processing or sent as fresh tissues as soon as possible. Long transport times degrade fresh tissues	<b>Tissues can be maintained 'fresh' in RS-I over ice @ 0 – 4 °C or room temperature (RT) for 12-72 hours before processing – fresh, FFPE and SF samples available.</b>
<b>Processing</b>	Immediate use of fresh tissues otherwise <u>LIMITED</u> as <b>SF</b> or <b>FFPE</b> samples	<b>UNLIMITED choice of techniques – fresh tissues, FFPE and SF are available for analysis</b>
<b>Analyses</b>	<u>LIMITED</u> by previous <b>SF</b> or <b>FFPE</b> procedures	<b>UNLIMITED choice</b>
<b>Outcomes</b>	Query altered transformation of molecular profiles	<b>Retains morphological and molecular conformation including functional activity</b>
<b>Prognostics</b>	<b><i>Inconsistent</i></b>	<b><i>Definitive</i></b>