

Title	Nucleic Acid Extraction from Tissue
SOP Code	SOP119_02
Effective Date	04-Jan-2016

Site Approvals

Name and Title (typed or printed)	Signature	Date dd/Mon/yyyy

1.0 PURPOSE

This Standard Operating Procedure (SOP) outlines standardized procedures for biorepositories to follow when extracting nucleic acids from tissue specimens. This SOP does not describe detailed safety procedures for handling Human Biological Materials (HBMs) or hazardous chemicals.

2.0 SCOPE

This is the suggested process for RNA extraction from snap frozen tissue; including, tissue frozen in Optimal Cutting Temperature (OCT) medium, and DNA extraction from snap frozen tissue.

3.0 RESPONSIBILITIES

This procedure applies to all biorepository personnel responsible for extracting RNA or DNA from tissue.

4.0 DEFINITIONS

See Glossary of Terms.

5.0 PROCEDURE

5.1 General Extraction Considerations

- 5.1.1 Due to the sensitivity of nucleic acid amplification technologies, precautions should be taken to avoid cross contamination of samples.
- 5.1.2 Avoid moistening the rim of the spin columns with pipette tips, and avoid touching the column with the pipette tip.
- 5.1.3 Always use aerosol-barrier tips.
- 5.1.4 Avoid cross-contamination after each vortex step, and briefly centrifuge the tubes to remove droplets that may be on the lids of the tubes.
- 5.1.5 Close the lids of the spin columns before placing in the micro-centrifuge.
- 5.1.6 Flow-through generated after each centrifugation step may contain hazardous materials, and must be disposed of appropriately.
- 5.1.7 Open only one spin column at a time, and avoid creating aerosols.
- 5.1.8 Do not use any plastic-ware or glassware without first eliminating RNase or DNase contamination.
- 5.1.9 Take care not to introduce RNase or DNase into the sample during or after the purification procedure.
- 5.1.10 It is optimal to use sterile RNase-free or DNase-free disposable vessels, and solutions while working with nucleic acids. Microbiological aseptic technique is always optimal when working with nucleic acids.
- 5.1.11 Wear latex or vinyl gloves while handling reagents, tubes, and samples to prevent RNase and DNase contamination from the skin or surface of the laboratory. Change gloves frequently.
- 5.1.12 Keeps tubes closed whenever possible.
- 5.1.13 Keep purified RNA on ice.
- 5.1.14 Keep samples frozen at or below -80° C for long term storage.

5.2 Extraction of RNA from Frozen Tissue

Note: Volumes indicated are recommendations only, and should be scaled according to the size of the tissue specimen. A commercial RNA extraction kit is suggested for RNA extraction from specimens that are small in size.

- 5.2.1 Have materials and equipment ready before starting the procedure. Have as many tubes, and cryovials as needed labelled and ready.
- 5.2.2 Homogenization: Keep tissues frozen at -80°C until homogenization.
- 5.2.3 Homogenize tissue samples in 1-mL of TRIZOL reagent per 50-100-mg of tissue using a glass-Teflon or power homogenizer. Alternate RNase free methods for homogenizing frozen tissue, can be used if a homogenizer is not available.
- 5.2.4 The sample volume should not exceed 10% of the volume of the TRIZOL used for homogenization.
- 5.2.5 Phase Separation: Incubate the homogenized samples for 5 minutes at room temperature, to permit complete dissociation of nucleoprotein complexes.
- 5.2.6 Add 0.2-mL of chloroform per 1-mL of TRIZOL reagent. Cap tubes securely, and shake tubes vigorously by hand for 15 seconds.
- 5.2.7 Incubate at room temperature for 2-3 minutes.
- 5.2.8 Centrifuge the sample at no more than 12,000 x g for 10 minutes in a refrigerated centrifuge (2-8°C).
- 5.2.9 Following centrifugation, the mixture separates into distinct phases: a lower red phenol/chloroform phase, an interphase, and a colourless upper aqueous phase.
- 5.2.10 RNA remains exclusively in the aqueous phase, and this phase is about 60% of the volume of the TRIZOL reagent used for homogenization.
- 5.2.11 RNA precipitation: Transfer the aqueous phase to a clean tube (the organic phase may be saved if isolation of DNA or protein is required from this same sample).
- 5.2.12 Precipitate the RNA from the aqueous phase by mixing with isopropyl alcohol. Use 0.5-mL of isopropyl alcohol per 1-mL of TRIZOL reagent used for the original homogenization.
- 5.2.13 Incubate at room temperature for 10 minutes.

- 5.2.14 Centrifuge the sample at 18,000 x g for 10 minutes in a refrigerated centrifuge (2-8°C).
- 5.2.15 The precipitated RNA forms a gel-like pellet on the lower side, and the bottom of the centrifuge tube.
- 5.2.16 Following centrifugation, remove the supernatant. .
- 5.2.17 Wash the RNA pellet once with 75% ethanol adding at least 1-mL of 75% ethanol per 1-mL of TRIZOL reagent that was initially used for homogenization.
- 5.2.18 Mix the sample by gentle vortex.
- 5.2.19 Centrifuge at 7500 x g for 5 minutes, in a refrigerated centrifuge (2-8°C).
- 5.2.20 Remove the supernatant and re-dissolve the precipitated RNA: Briefly dry the RNA pellet. Air dry or vacuum dry for 5-10 minutes. Take care not to completely dry the pellet, as this will cause difficulty with dissolving the RNA.
- 5.2.21 Dissolve the pellet in an appropriate volume of RNase free water.
- 5.2.22 Store the dissolved RNA at -80°C or lower.
- 5.2.23 Record the storage location.

5.3 Extraction of RNA from Tissue Frozen in OCT

Note: Volumes indicated are recommendations only, and should be scaled according to the size of the tissue sample. A commercial RNA extraction kit is suggested for RNA extraction from OCT sections.

- 5.3.1 Have materials and equipment ready before starting the procedure. Have as many tubes, and cryovials as needed labelled and ready.
- 5.3.2 Take several (5-10) 3-µm OCT sections using a cryostat, and place them in a pre-cooled microfuge tube. Ensure that the sections do not thaw before the next step.
- 5.3.3 Follow the detailed procedure outlined in the RNA extraction kit manual.
- 5.3.4 Immediately after the procedure, place extracted and re-suspended RNA on ice.
- 5.3.5 Store extracted RNA at or below -80°C.

5.3.6 Record the storage location.

5.4 Extraction of DNA from Frozen Tissue

5.4.1 Have materials and equipment ready before starting the procedure. Have as many tubes, and cryovials as needed labelled and ready.

5.4.2 Cut tissue (mince) into small pieces with sterile scissors or scalpel blade. Alternatively, wrap frozen tissue in aluminum foil, and fragment with a hammer.

5.4.3 Follow the detailed procedure outlined in the DNA extraction kit manual.

5.4.4 Store DNA at 4°C, or for long term storage at -20° or below.

6.0 REFERENCES

Health Canada, Food and Drug Regulations, Part C, Division 5, Drugs for Clinical Trials Involving Human Subjects, (Schedule 1024), June 20, 2001.

Health Canada, Guidance for Industry, Good Clinical Practice: Consolidated Guideline, ICH Topic E6, 1997.

2011 NCI Best Practices for Specimen Resources. Office of Biorepositories and Biospecimen Research, National Cancer Institute, Bethesda, MD.
<http://biospecimens.cancer.gov/bestpractices/2011-NCIBestPractices.pdf>

ISBER Best Practices for repositories: Collection, storage, retrieval and distribution of biological materials for research, 3rd Edition, <http://www.isber.org>

CTRNET Standard Operating Procedures, Canadian Tissue Repository Network

7.0 REVISION HISTORY

SOP Code	Effective Date	Summary of Changes
SOP119_01	01-Aug-2012	Original version
SOP119_02	04-Jan-2016	2.0 Scope, 5.1.4, 5.1.6, 5.1.10, 5.1.11, 5.1.14, 5.2.3, 5.2.6, 5.2.9, 5.2.12, 5.3.5, 5.4.4: rewording for clarification 5.2.14: correction in centrifuge speed. 5.2.15 – 5.2.23: rewording/renumbering for clarification. Updated references. Removed OTRN logo.