


Saliva Collection:



1. Open the package containing the saliva collection device. Place contents on a clean and dry surface.
 2. Place lips over the opening of the collection Tube Funnel. Collect saliva until the sample level reaches the red indicator line. (Ensure that the saliva reaches the red indicator line, not just the foam.)
 3. When collection has been completed, unscrew the Funnel from the collection tube by rotating the tube in the direction shown until the tube is completely detached from the funnel.
 4. Remove the cap from the bottom of the collection tube and screw it tightly closed.
 5. Gently invert the vial 20 times to stabilize the sample.
 6. Remove the barcode label from the outside of the saliva device package and apply to the collection tube.
 7. Place collection tube into biohazard bag with absorbent pad. Ensure the seal is closed completely. *Note: The thin "paper" sheet in the biohazard bag is the absorbent pad. Please leave in biohazard bag.*
- ✓ After specimens are collected, please place all biohazard bags with specimens into the kit box.
 - ✓ Please ensure the barcode label has been applied to the Saliva Collection Tube.
 - ✓ Please ensure your Date of Birth has been written on each barcode label in the "DOB: _____" location.

Leftover supplies can also go into the box or be discarded in the trash.